



RetenSpect²⁵



Focused on force, flexibility and retention


➔ Key Features

Unparallel load verification and superior stiffness and alignment



Single column tabletop system for force testing

Graphical representation of load, time, motion and optional pictorial data



High capacity motor for better reliability

➔ Specifications

- ➔ Force measurement capacity: 25 N (other options available)
- ➔ Force measurement resolution: 0.01 N
- ➔ Z Axis traveling capacity: 180 mm
- ➔ Z Axis traveling resolution: 1μ

Distributed by:

