



# RetenSpect<sup>25</sup>



**Focused on force, flexibility and retention**


## ➔ Key Features

**U**nparallel load verification and superior stiffness and alignment



**S**ingle column tabletop system for force testing

**G**raphical representation of load, time, motion and optional pictorial data




**H**igh capacity motor for better reliability

## ➔ Specifications

- ➔ Force measurement capacity: 25 N (other options available)
- ➔ Force measurement resolution: 0.01 N
- ➔ Z Axis traveling capacity: 180 mm
- ➔ Z Axis traveling resolution: 1μ

Distributed by:



**SAHAJANAND**  
TECHNOLOGIES (P) LTD.

Where **TECHNOLOGY**  
Becomes **REVOLUTION**

"Sahajanand House", Parsi Street, Saiyedpura, Surat – 395 003

Tel : +91 – 261 – 2451451 – 58 • Fax : +91 – 261 – 2427947

Web : [www.stentspect.com](http://www.stentspect.com)

E-mail : [stentinfo@stentspect.com](mailto:stentinfo@stentspect.com)